RZera User Manual

CSNS Data Analysis Group

RZera User Manual

The purpose of this manual is to explain to users the workflow for data reduction using RZera, and to guide users step by step through the data reduction process. If the data in your paper have been processed using RZera, please cite reference [1]. If you use and adopt the prediction results of the AI module "AIndex" in RZera in your paper, please cite reference [2]. If you use the small-angle scattering module, please cite reference [3]. If you encounter any issues during the use of the diffraction reduction function, please contact wanghao96@ihep.ac.cn. For issues during the use of the small-angle reduction function, please contact zhongjiajun@ihep.ac.cn. For collaboration inquiries, please contact durong@ihep.ac.cn.

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1. Instructions for Using Powder Diffraction Spectrometer Software

1.1 Data Processing Workflow for Powder Diffraction in RZera

RZera splits the data reduction workflow for neutron powder diffraction into several sequential steps, as shown in Figure 1.1-1. In the first step, the data of the sample, sample cell, V/VNi, V/VNi cell, and empty background are each subjected to time focusing, which converts the two-dimensional histogram data into one-dimensional I-d curves. RZera provides a dedicated "Time Focusing" page for this operation. In the second step, the one-dimensional curves are used for background subtraction and vanadium normalization to obtain the reduced I-d curve of the sample. This operation is handled by the "Reduction" page in RZera. Note that some spectrometers do not require empty background data; in such cases, the data processing steps related to the empty background can be omitted. In the subsequent figures, the background refers specifically to the cell data. After that, users can continue to process the data to obtain the PDF. In RZera, there are two methods for obtaining the PDF: one is the Merge + Fourier transform method, and the other is the Hermite polynomial fitting method.


Figure 1.1-1 Data reduction workflow in RZera without empty background subtraction. In the figure, "background" refers to the cell. If empty background subtraction is required, the sample, V/VNi, and cell data must first be subtracted once by the empty background data.

In the above reduction workflow, almost every step supports visualization, allowing users to conveniently check the changes in the data throughout the entire process. However, because the steps are divided in considerable detail, the operation can be rather tedious. To address this, RZera also provides a quick reduction function: once the relevant parameters are properly configured, batch data reduction can be achieved with a single click, as shown in Figure 1.1-2. First, users still need to perform time focusing on the V/VNi data and its background data, but these two datasets can be reused across a beamtime experiment. Then, RZera provides a "Quick Reduction" page, which accepts the time-focused V/VNi and V/VNi background, as well as the pre-time-focusing sample and sample background (post-time-focusing background data are also acceptable) as inputs. After the user has set all parameters on the interface, the reduced sample data can be obtained with one click. RZera also offers a configuration saving feature: when data from different samples are collected next time, users only need to import the saved configuration, replace the sample data, and then perform the one-click

reduction again.


Figure 1.1-2 Quick data reduction workflow in RZera.

1.2 Data Reduction Using Histogram Data

1.2.1 File Preparation

Users need to prepare four files, namely the nxs files (provided by the spectrometer) for the V or VNi data, the V or VNi background, the sample data, and the sample background. The V or VNi background and the sample background may be the same dataset; in that case, only three files need to be prepared. As shown in Figures 1.2-1 and 1.2-2, there are three folders named after RUN numbers, which contain the nxs files for the background, V, and Si standard sample, respectively.

RUN0023566	→ the background of Si sample and V	2025/7/3 16:26	文件夹
RUN0023567	→ V	2025/5/21 13:48	文件夹
RUN0023568	→ Si Sample	2025/7/3 16:26	文件夹

Figure 1.2-1 Folders storing nxs files, with folder names corresponding to RUN numbers. The RUN number is used to label each data collection experiment.

Figure 2 Data inside a folder named by the RUN number. The nxs file is "detector.nxs" as shown in the figure, which stores the histogram data and some metadata of the experiment. There may be other data files in this folder as well, but this chapter, "Data Reduction Using Histogram Data," only uses the "detector.nxs" file.

1.2.2 Entering the RZera data reduction interface for the corresponding spectrometer


Figure 1.2-3 How to enter the data reduction window for a specified spectrometer.

As shown in Figure 1.2-3, after opening RZera, there is a drop-down box for selecting a spectrometer. Once the spectrometer used in the experiment is selected, the corresponding data reduction window for that spectrometer will pop up.

1.2.3 Performing Time Focusing on the V/VNi Data

The screenshot shows the 'Load and Time Focusing' window with the following fields and steps indicated by red arrows:

- 1. Select Histogram Data:** Points to the 'Histogram Data' tab.
- 2. Click Select Run, select target nxs file:** Points to the 'Select Run' button and the text 'RUN0023567/detector.nxs'.
- 3. Click Select Detector, select target detector:** Points to the 'Select Detector' dropdown menu showing 'group2; group3; group4; group5; group6; group7'.
- 4. Click Instru Folder, select folder with instrument files in it:** Points to the 'Instru Folder' text field containing ':search_or_work/qt_develop/rzera_offline/param_data/BL16/pidInfo'.
- 5. Check Crop Data, change wave min and wave max:** Points to the 'Crop Data' checkbox (checked) and the 'Wave Min (Å)' field (0.1) and 'Wave Max (Å)' field (4.5).
- 5. Check Carpenter Correction, make the absorption and multiple scattering correction, fill the below experimental information:** Points to the 'Carpenter Correction' checkbox (checked).
- 6. Check offset, then click "Offset Folder", select folder with offset files in it:** Points to the 'Offset' checkbox (checked) and the 'Offset Folder' text field containing ':search_or_work/qt_develop/rzera_offline/param_data/BL16/offsetfiles_RUN0023568_RT_HOT04'.

Other visible fields include 'Mask Data' (unchecked), 'Mask Folder', 'Sample Name' (V:1), 'Example: Si:1-O:2', 'Sample Mass (g)' (11.166), 'Sample Radius (cm)' (0.45), 'Beam Height (cm)' (3.0), 'Sample Height (cm)' (3.0), 'Calculated Num Density: 0.069140 Å⁻³', and 'Inputed Num Density (Å⁻³)(priority)' (0.0721). Buttons at the bottom include 'Merge Load', 'Batch Load', and 'Delete'.

Figure 1.2-4 Information input steps in "Load and Time Focusing"

In the pop-up reduction window, there is a Time Focusing page, where the first module to be operated is "Load and Time Focusing", as shown in Figure 4. When importing data following the instructions in Figure 1.2-4, several points should be noted:

(1) Regarding the "Batch Load" option: Users can import data from multiple Runs simultaneously. If "Batch Load" is checked, RZera will perform time focusing on each of the imported datasets separately. However, if it is not checked, RZera will merge the data from these Runs and then perform time focusing on the merged data.

(2) Regarding "Instru Folder", RZera will provide this.

(3) Regarding "Sample_Name" under "Carpenter Correction": Enter the ratio according to the actual molar fraction of the V:Ni sample. If the exact number density value is known, instead of using the calculated value, you may check "Inputed Num Density" and manually enter the accurate value.

(4) Regarding "Offset Folder", RZera or the spectrometer team can provide this.

After completing the information input, click "Load", and the software will begin loading the data and performing time focusing based on the entered information. The loaded data will then be saved in the list at the bottom of this module, as shown in Figure 1.2-5.

Merge Load Batch Load Delete			
	Name	Detector	Select
1	RUN0023567_crop_carpenterCorrection_offset_time(2025-05-28-15:59:16)	group2	☐
2	RUN0023567_crop_carpenterCorrection_offset_time(2025-05-28-15:59:20)	group3	☐
3	RUN0023567_crop_carpenterCorrection_offset_time(2025-05-28-15:59:24)	group4	☐
4	RUN0023567_crop_carpenterCorrection_offset_time(2025-05-28-15:59:31)	group5	☐

Figure 1.2-5 Data obtained after loading.

After the data are loaded, this page also provides peak removal and smoothing functional modules, located below the "Load and Time Focusing" module, as shown in Figure 1.2-6. These modules are collapsed by default and are used to remove unwanted diffraction peaks from the data. For vanadium, we need to remove all of its diffraction peaks.

	Name	Detector	Select
3	RUN0023567_crop_carpenterCorrection_offset_time(2025-05-28-15:59:24)	group4	☐
4	RUN0023567_crop_carpenterCorrection_offset_time(2025-05-28-15:59:31)	group5	☐
5	RUN0023567_crop_carpenterCorrection_offset_time(2025-05-28-15:59:35)	group6	☐
6	RUN0023567_crop_carpenterCorrection_offset_time(2025-05-28-15:59:39)	group7	☐

► Cut Peaks

► Smooth

► Save Data

Figure 1.2-6 Peak removal and smoothing functional modules.

After clicking "Cut Peaks", the module will expand, as shown in Figure 1.2-7. This module consists of two parts: the upper part is used to automatically calculate the positions of diffraction peaks to be removed based on a CIF file, and the lower part is for manually locating diffraction peak positions.

Cut peaks with cif

▼ Cut Peaks
☐ Plot

Add Cif
Delete
Check
Add Peaks

Name	a	b	c	α	β	γ	A	B

Add Default Peak
Add Peak
Delete

Detector	Peak Left	Peak Right	Select
Cut peaks manually			

Figure 1.2-7 The "Cut Peaks" functional module.

RZera has built-in default peak positions for vanadium. Simply click "Add Default Peak" in the manual peak removal section to automatically import the V peak positions, as shown in Figure 1.2-8.


Figure 1.2-8 Importing default V peaks. (a) After clicking the button, select the type of default peaks to import. (b) Select the detector for which the peak information is to be imported. The detectors in this list are consistent with those selected during the previous data loading step. (c) The peak position information displayed after import.

After importing, check the "plot" boxes in both the "Load and Time Focusing" module and the "Cut Peaks" module, then click the "Run & Plot" button at the bottom of the entire page. The "Run & Plot" button will execute the enabled functional modules in order from top to bottom. Additionally, if the "plot" box in a module is checked, the data processed by that module will be added to the plotting window and the plotting

window will pop up, enabling data visualization, as shown in Figure 1.2-9.


Figure 1.2-9 Comparison of data before and after peak removal. The reason for this comparison is that the software executes modules from top to bottom: first, "Load and Time Focusing" is executed, and since the "plot" box in this module is checked, the data before peak removal are saved in the plotting window. Then, "Cut Peaks" is executed, and after peak removal, since the "plot" box in this module is also checked, the data after peak removal are saved in the plotting window.

After peak removal, the removed portions are filled with straight lines. RZera provides the "Smooth" module (Figure 1.2-6) to smooth the data, and its expanded view is shown in Figure 1.2-10. Clicking "Add Default Smooth" within this module allows you to apply the built-in smoothing parameters of RZera with one click.


Figure 1.2-10 (a) The expanded interface of the "Smooth" module. (b) The interface after clicking "Add Default Smooth".

After checking the "plot" box in the "Smooth" module as well, click the "Run & Plot" button to achieve a comparison of the data before and after smoothing, as shown in Figure 1.2-11.


Figure 1.2-11 Comparison of data from detector group 7 before and after smoothing. It can be seen that the data are significantly smoother after smoothing.

After smoothing, expand the final "Save Data" module (as shown in Figure 1.2-12) to prepare for saving the processed data. After configuring the "Save Data" module, click "Run & Plot" again, and the processed data will be saved to the specified path.


Figure 1.2-12 The interface of the "Save Data" module. The "Data Type" must be consistent with the type of data currently being processed. In this example, we are processing V or VNi sample data, so "V" should be selected. Then click "Save Path" to specify the save location.

The processed file is a ".nc" file, which stores the processed data and the processing history internally, as shown in Figure 1.2-13.

	v_RUN0023567_crop_carpeneterCorrection_offset_time(2025-05-28-173126)_group2.nc	2025/5/28 17:33
	v_RUN0023567_crop_carpeneterCorrection_offset_time(2025-05-28-173130)_group3.nc	2025/5/28 17:33
	v_RUN0023567_crop_carpeneterCorrection_offset_time(2025-05-28-173135)_group4.nc	2025/5/28 17:33
	v_RUN0023567_crop_carpeneterCorrection_offset_time(2025-05-28-173146)_group5.nc	2025/5/28 17:33
	v_RUN0023567_crop_carpeneterCorrection_offset_time(2025-05-28-173151)_group6.nc	2025/5/28 17:33
	v_RUN0023567_crop_carpeneterCorrection_offset_time(2025-05-28-173157)_group7.nc	2025/5/28 17:33

Figure 1.2-13 Saved Data

1.2.4 Performing Time Focusing on the V/VNi Background

The operation procedure for this step is the same as in the previous step (Step 3), but there are several differences or points to note, as follows:

- (1) The nxs file to be loaded should be replaced with the background


corresponding to V/VNi.

(2) If Carpenter Correction was applied to the V/VNi data during "Load and Time Focusing", it is recommended to apply exactly the same Carpenter Correction to the background data as well, in order to cancel out its effect during the subsequent background subtraction.

(3) Since the background may not contain diffraction peaks, the Cut Peaks and Smooth steps may not be necessary, depending on the actual situation.

(4) When saving the data at the end, "Data Type" should be set to "V Bkg".

The final saved data are shown in Figure 1.2-14.

vBG_RUN0023566_crop_carpenterCorrection_offset_time(2025-06-03-142415)_group2.nc
vBG_RUN0023566_crop_carpenterCorrection_offset_time(2025-06-03-142418)_group3.nc
vBG_RUN0023566_crop_carpenterCorrection_offset_time(2025-06-03-142423)_group4.nc
vBG_RUN0023566_crop_carpenterCorrection_offset_time(2025-06-03-142431)_group5.nc
vBG_RUN0023566_crop_carpenterCorrection_offset_time(2025-06-03-142434)_group6.nc
vBG_RUN0023566_crop_carpenterCorrection_offset_time(2025-06-03-142439)_group7.nc

Figure 1.2-14 The saved V/VNi background data.

1.2.5 Performing Time Focusing on the Sample

The operation procedure for this step is the same as Step 3, but there are several differences or points to note, as follows:


(1) The nxs file to be loaded should be replaced with that of the sample.

(2) Carpenter Correction is not necessarily required for the sample. If only diffraction data are needed, it can be skipped. However, if further processing to obtain the PDF is required, then it is needed. If the correction is to be applied, it is necessary to ensure that the sample is either a cylinder or a powder contained in a cylindrical container, and that all relevant input parameters are entered correctly; otherwise, the correction may be inaccurate.

(3) Since the sample may not contain extraneous diffraction peaks, the Cut Peaks and Smooth steps may not be needed, depending on the specific situation.

(4) When saving the data at the end, "Data Type" should be set to "Sample".

The final saved data are shown in Figure 1.2-15.


```

sam_RUN0023568_crop_carpenterCorrection_offset_time(2025-06-03-143746)_group2.nc
sam_RUN0023568_crop_carpenterCorrection_offset_time(2025-06-03-143749)_group3.nc
sam_RUN0023568_crop_carpenterCorrection_offset_time(2025-06-03-143754)_group4.nc
sam_RUN0023568_crop_carpenterCorrection_offset_time(2025-06-03-143801)_group5.nc
sam_RUN0023568_crop_carpenterCorrection_offset_time(2025-06-03-143805)_group6.nc
sam_RUN0023568_crop_carpenterCorrection_offset_time(2025-06-03-143809)_group7.nc

```

Figure 1.2-15 The saved sample data.

1.2.6 Performing Time Focusing on the Sample Background

The operation procedure for this step is the same as Step 3, but there are several differences or points to note, as follows:

- (1) The nxs file to be loaded should be replaced with that of the sample background.
 - (2) If Carpenter Correction was applied to the sample during "Load and Time Focusing", it is recommended to apply exactly the same Carpenter Correction to the background data as well, in order to cancel out its effect during the subsequent background subtraction.
 - (3) Since the background may not contain diffraction peaks, the Cut Peaks and Smooth steps may not be necessary, depending on the actual situation.
 - (4) When saving the data at the end, "Data Type" should be set to "Sample Bkg".
- The final saved data are shown in Figure 1.2-16.

```

samBG_RUN0023566_crop_carpenterCorrection_offset_time(2025-06-03-151459)_group2.nc
samBG_RUN0023566_crop_carpenterCorrection_offset_time(2025-06-03-151503)_group3.nc
samBG_RUN0023566_crop_carpenterCorrection_offset_time(2025-06-03-151507)_group4.nc
samBG_RUN0023566_crop_carpenterCorrection_offset_time(2025-06-03-151515)_group5.nc
samBG_RUN0023566_crop_carpenterCorrection_offset_time(2025-06-03-151518)_group6.nc
samBG_RUN0023566_crop_carpenterCorrection_offset_time(2025-06-03-151522)_group7.nc

```

Figure 1.2-16 The saved sample background data.

1.2.7 Using the data obtained in Sections 1.2.3–1.2.6 to obtain the reduced diffraction pattern.

After completing Steps 3–6, we have obtained the time-focused data for four types: sample, sample background, V/VNi, and V/VNi background. Next, we need to switch

the tab page from "Time Focusing" to "Reduction". The interface of the Reduction page is shown in Figure 1.2-17.


Figure 1.2-17 The "Reduction" tab page.

First, we operate the "Load Focused Data" module. Click "Load" to import all the data obtained in Steps 3–6, as shown in Figure 1.2-18.


Figure 1.2-18 (a) Schematic diagram of the "Load Focused Data" module operation. (b) The "Load Focused Data" module after data loading.

Next, consider the "Normalization" module, which offers two options: "By Proton Charge" and "By Monitor Counts". These refer to normalizing the data by proton charge and by monitor counts, respectively. In theory, either choice is correct; however, there may be cases where a monitor was not used during the experiment, in which case selecting "By Monitor Counts" will cause the program to report an error. Therefore, the appropriate option should be chosen based on the actual experimental conditions.

Next, expand the "d_Rebin" module and click "Add Default Rebin". After selecting the detectors to be processed in the pop-up window, the module will load the

default d_rebin parameters for each detector, as shown in Figure 1.2-19. After loading the default parameters, users can also modify the d_rebin parameters directly in the interface.


Figure 1.2-19 (a) The workflow for operating the "d_Rebin" module. (b) The "d_Rebin" module after loading the default parameters.

Next, operate the subsequent "Background Subtraction" module as shown in Figure 1.2-20. This module is dedicated to background subtraction. It can be noted that there are multiple consecutive "Background Subtraction" modules in the interface, which are used for subtracting the background from the sample and V/VNi data respectively (there is no required order for the background subtraction). Click "Load Sample" to load the sample or V/VNi data, and then click "Load Background" to load the corresponding background. After loading, the module provides real-time display of the data after background subtraction. The "Background Scale" in the module refers to the scale factor multiplied to the background before subtraction. Sliding the "Background Scale" slider allows adjustment of this factor, and the changes in the data after applying the factor are reflected in real time in the plot on the right. The default scale factor is 1, but there may be cases where the background is higher than the data; in such cases, the scale factor must be adjusted to ensure that all data remain positive after background subtraction. For guidance on which value to choose as the scale factor, please consult the spectrometer scientist for the corresponding instrument.


Figure 1.2-20 The interface after loading the sample, V/VNi, and their corresponding backgrounds. After the data are loaded, the right panel will display the background-subtracted data. If the data do not appear, you can click the "home" button (the small black house icon) in the lower-left corner of the plotting area, and the data will be displayed normally.

Next, operate the subsequent "Division" module. Click "Load Sample" to import the sample data, and click "Load Calibration" to import the V/VNi data. The data imported here will be divided after having undergone background subtraction from the preceding modules. The module after data import is shown in Figure 1.2-21.

Division		☐ Plot
Load Sample	5-143805_group6; sam_RUN0023568_crop_carpen...Correction_offset_time(2025-06-03-143809)_group7	
Load Calibration	28-173151_group6; v_RUN0023567_crop_carpen...Correction_offset_time(2025-05-28-173157)_group7	

Figure 1.2-21 The "Division" module after importing the sample and V/VNi data.

Next is a skippable module, "AIndex", which is developed based on research results published in 2025 by the China Spallation Neutron Source data analysis group. Its function is AI-based prediction of space groups and unit cell parameters. Users only need to import the reduced data to predict the space group and unit cell parameters with one click, and can easily check whether the prediction results are reliable.

First, click the "Load Pattern" button to import the reduced sample data, as shown in Figure 1.2-22 (a)(b). Then, click the "Prediction" button, and the software will call the trained deep neural network to provide prediction results for the crystal system,

space group, and unit cell parameters corresponding to each diffraction pattern, which are displayed in the list panel on the right, as shown in Figure 1.2-23. The panel contains two columns: the first column shows the prediction results, and the second column shows the confidence level of each prediction. By selecting a prediction result in the list panel and clicking the "Check" button, AIndex will calculate the diffraction peak positions using the predicted results and display them together with the diffraction pattern, as shown in Figure 1.2-23, allowing users to conveniently verify the accuracy of the prediction.

There are certain limitations to using this module, as well as specific requirements for the data: (1) Currently, AIndex can only perform predictions for cubic, tetragonal, trigonal, and hexagonal crystal systems. (2) The data must be from a single-phase sample; otherwise, correct prediction cannot be performed. (3) The unit cell parameters of the sample should be within the range of 3–15 Å; otherwise, correct prediction cannot be performed. (4) The data should ideally cover the range 0.5–10 Å, or it must be ensured that no diffraction peaks exist in the uncovered regions; otherwise, the prediction results are likely to be inaccurate.



Figure 1.2-22 (a) The workflow for importing data into "AIndex". (b) After data import, AIndex reorganizes the data into a diffraction pattern consisting of 4196 equally spaced points over the range 0.5–10 Å and displays it. This reorganization is intended to match the input format required by the neural network and is used only internally within this module; it does not affect the final reduction results.

1. Click "AI Prediction".

Figure 1.2-23 The workflow for predicting space group and unit cell parameters using AIndex with the assistance of AI.

Then, operate the "Save Reduced Data" module. This module is very simple: after expanding it, just click "Save Path" to specify the save location for the data. Finally, click "Run & Plot", and the software will process the modules one by one from top to bottom. The processed data will then be saved to the path specified in the "Save Reduced Data" module, as shown in Figure 1.2-24.


Figure 1.2-24 The saved reduced data. For each dataset, four files are saved: time-of-flight data files for data processing in the three software packages "FullProf", "GSAS", and "Z-Retiveld", as well as a txt file containing the d-spacing data.

1.2.8 Quick Reduction

After obtaining the time-focused data of V/VNi and its cell (container) through Sections 1.2.3 and 1.2.4, quick data reduction can be performed on the Quick Reduction page in RZera. The following figure shows the interface of Quick Reduction along with configuration and usage instructions.

The screenshot shows the 'MPI' Quick Reduction interface. At the top, there's a 'MainWindow' title bar and a 'Reload' button. Below it, a status bar says 'Instrument configure file has been loaded. Beamline name BL16:MPI'. The main section is titled 'Quick Reduction' and contains several tabs: 'Histogram Sample', 'Event Sample', 'Select Run', 'Batch Load', 'Select Sample Cell', 'Select Van', 'Select Van Cell', 'Select Detector', 'Instru Folder', 'Mask Data', 'Crop Data', 'Sample Correction', 'Offset', 'Normalization Methods', 'Bkg Subtraction', and 'Plot'. Red arrows point to various fields with instructions: 'select the sample runs (nxs)' points to 'Select Run'; 'if check, deal the runs seperately, otherwise sum them before dealing' points to 'Batch Load'; 'sample cell data, nxs/nc optional' points to 'Select Sample Cell'; 'V/VNi nc data after time focusing' points to 'Select Van'; 'V/VNi cell nc data after time focusing' points to 'Select Van Cell'; 'target detectors' points to 'Select Detector'; 'folder with instrument files in it' points to 'Instru Folder'; 'check it for wave truncation' points to 'Wave Min (Å) 0.1'; 'folder with offset files in it' points to 'Offset Folder'; 'Normalization Methods' points to 'By Proton Charge' and 'By Monitor Counts'; 'Scale factor of sample and V/VNi cell' points to 'Sample Bkg Scale 1.0' and 'Van Bkg Scale 0.5'; and 'click the "reduction" finally' points to the 'Reduction' button. The bottom of the interface has three buttons: 'Run & Plot', 'Save Configure', and 'Load Configure'.

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The screenshot shows a software interface for adjusting parameters. It is divided into three main sections, each with a title and several input fields:

- Transmittance**
 - Brag threshold: 0.04
 - Min threshold: 0.7
- Orientation**
 - Residual standard: 0.4
 - Count threshold: 5
 - Angle range: 30.0
 - Q grid: 82
 - Radius: 10
- Wave slice**
 - Deviation threshold: 0.2
 - Accept ratio: 0.8

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Output directory		
	Configuration_initial.json	} Input Files
	Sample_info.xlsx	
	1D IQ_runno_SubEC/EB_PC/M2_batch.dat	Sample Iq after correction and bkg subtraction
	1DRaw IQ_runno_Raw_PC/M2_batch.dat	Sample and bkg scattering curve after correction
	2D IQ2D_runno_SubEC/EB_PC/M2_batch.dat	Sample 2D data after correction and bkg subtraction
	2DRaw IQ2D_runno_Raw_PC/M2_batch.dat	Sample and bkg 2D scattering curve after correction
	Trans Trans_all/1D_PC/M2_runno_batch.dat	Transmission curve
	Check IQ_runno_check_no.dat	Iq with multi-wave check
	DetectorImage RUN00runno_detector_image.png	2D detector picture before correction
	logs sans_reduction_date.log	log
		runno_trans_check.png
	auto_reduction Quality_runno_SubEC/EB_PC/M2_batch.dat	the result of auto check for reduction